

Product range: Serologic Water Bath

Polycarbonate Water Bath – 3 litre Capacity

Reference: BMP30

General Specifications:

Net capacity: 3 litres

Temperature range: ambient +7°C to +60°C

Construction:

Polycarbonate frame

External dimensions W x H x D: 360 x 120 x 130 mm

Internal dimensions W x H x D: 280 x 110 x 80 mm

Gross weight: 3kg

Net weight: 1,5kg

Factory setpoint: +37°C and +56°C

Water lack safety



Options:

Polycarbonate lid

1 rack for 60 tubes diam. 14mm

Thermometer

Electric data:

Power supply: Voltage 230V ~ +/-10%, 50 Hz

Energy consumption: 250 W

Warranty: 2 years parts

Cleaning and decontamination with detergents and antimicrobials used in hospitals
(Ammonium quaternary / except pure bleach)

Device complies with CE specifications

Designed in France (Meyzieu 69) and assembled in Romania