

Product range: Serologic Water Bath

Polycarbonate Water Bath – 6 litre Capacity

Reference: BMP60

General Specifications:

Net capacity: 6 litres

Temperature range: ambient +7°C à +60°C

Construction:

Polycarbonate frame

External dimensions W x H x D: 370 x 150 x 210 mm

Internal dimensions W x H x D: 290 x 110 x 190 mm

Gross weight: 4kg

Net weight: 2kg

Factory setpoint: +37°C and +56°C

Water lack safety



Options:

Polycarbonate lid

2 racks for 60 tubes diam. 14mm

Thermometer

Electric data:

Power supply: Voltage 230V ~ +/-10%, 50 Hz

Energy consumption: 250 W

Warranty: 2 years parts

Cleaning and decontamination with detergents and antimicrobials used in hospitals
(Ammonium quaternary / except pure bleach)

Device complies with CE specifications

Designed in France (Meyzieu 69) and assembled in Romania