






Deep Freezer's storage elements

Side access aluminium racks	
	BM/EL3E690AC : 3-level aluminum rack with dividers
	External dimensions l x p x h : 140 x 730 x 343 mm Useful height of each level : 10,9 cm Capacity in boxes: 15
	BM/EL4E690AC : 4-level aluminum rack with dividers
	External dimensions l x p x h : 140 x 730 x 343 mm Useful height of each level : 8,0 cm Capacity in boxes: 20
	BM/EL5E690AC : 5-level aluminum rack with dividers
	External dimensions l x p x h : 140 x 730 x 343 mm Useful height of each level : 6,3 cm Capacity in boxes: 25
	BM/EL6E690AC : 6-level aluminum rack with dividers
	External dimensions l x p x h : 140 x 730 x 343 mm Useful height of each level : 5,4 cm Capacity in boxes: 30
	BM/TIR : stainless steel bulk storage drawer
	External dimensions l x p x h : 140 x 730 x 343 mm 4 drawers per compartment


Front access racks with drawers	
	BM/EL2T690A : 2-drawers aluminum front access rack
	External dimensions : l x p x h : 288 x 730 x 344 mm Inner dimension of each drawer l x p x h : 280 x 705 x 164 mm Useful height of each drawer : 16,4 cm / 2 drawers per rack
	BM/EL3T690A : 3-drawers aluminum front access rack
	External dimensions : l x p x h : 288 x 730 x 344 mm Inner dimension of each drawer l x p x h : 280 x 705 x 109 mm Useful height of each drawer : 10,9 cm / 3 drawers per rack
	BM/EL4T690A : 4-drawers aluminum front access rack
	External dimensions : l x p x h : 288 x 730 x 344 mm Inner dimension of each drawer l x p x h : 280 x 705 x 80 mm Useful height of each drawer : 8 cm / 4 drawers per rack
	BM/EL5T690A : 5-drawers aluminum front access rack
	External dimensions l x p x h : 288 x 730 x 344 mm Inner dimension of each drawer l x p x h : 280 x 705 x 63 mm Useful height of each drawer : 6,3 cm / 5 drawers per rack
	BM/EL6T690A : 6-drawers aluminum front access rack
	External dimensions l x p x h : 288 x 730 x 344 mm Inner dimension of each drawer l x p x h : 280 x 705 x 53 mm Useful height of each drawer : 5,3 cm / 6 drawers per rack

Rack capacity of Froilabo's Deep Freezers




340L

Rack capacity : 8




515L

Rack capacity : 12



690L

Rack capacity : 16



1000L

Rack capacity : 24

For the following racks :

BM/EL3E690AC


BM/EL4E690AC

BM/EL5E690AC

BM/EL6E690AC


BM/TIR

Rack with drawers capacity in Froilabo's Deep Freezers




340L

Rack capacity : 6




515L

Rack capacity : 6



690L

Rack capacity : 8



1000L

Rack capacity : 12

For the following racks with drawers :

BM/EL2TE690A

BM/EL3TE690A

BM/EL4TE690A

BM/EL5TE690A

BM/EL6TE690A