

# Thermostatic Water Bath

## BMP 60



### Specially designed for serology

Ideal for modern laboratories, our serological bath is available in 2 volumes; 3 and 6 litres, and consists of a transparent plexiglass and a control interface via a button. This makes it very easy to use and set the ambient temperature between 7°C to 60°C.

#### PERFORMANCE


Ambient temperature range	7°C to 60°C
Factory preset	37°C and 56°C
Capacity	6L

#### WARRANTY

Parts	2 years
-------	---------



Water lack  
safety



Available in  
2 volumes

# Technical Specifications

SEROLOGIC WATER BATH	
MODEL	BMP 60
<b>GENERAL SPECIFICATIONS</b>	
Capacity (net)	6 liters
Temperature Range (ambient)	7°C to 60°C
Construction	Polycarbonate frame
External Dimensions (W x H x D)	370 x 150 x 210 mm
Internal Dimensions (W x H x D)	290 x 110 x 190 mm
Weight (gross)	4kg
Weight (net)	2kg
Factory Setpoint	37°C and 56°C
Water lack safety	Yes
Warranty	2 years
<b>ELECTRICAL DATA</b>	
Power supply	230V ~ +/-10%, 50 Hz
Energy consumption	250 W
<b>COMPLIANCE</b>	
Device complies with CE specifications	
<b>OPTIONS</b>	
Polycarbonate lid	
2 racks for 60 tubes diam. 14mm	
Thermometer	

\* Cleaning and decontamination with detergents and antimicrobials used in hospitals (Ammonium quaternary / except pure bleach)